Software

Open Access

WASP: a Web-based Allele-Specific PCR assay designing tool for detecting SNPs and mutations

Pongsakorn Wangkumhang¹, Kridsadakorn Chaichoompu¹, Chumpol Ngamphiw¹, Uttapong Ruangrit¹, Juntima Chanprasert¹, Anunchai Assawamakin² and Sissades Tongsima^{*1}

Address: ¹Biostatistics and Informatics Laboratory, Genomics Institute, National Center for Genetic Engineering and Biotechnology, Thailand Science Park, Pathumtani, Thailand and ²Division of Molecular Genetics, Department of Research and Development, Faculty of Medicine Siriraj Hospital, Bangkok, Thailand

Email: Pongsakorn Wangkumhang - pongsakorn.wan@biotec.or.th; Kridsadakorn Chaichoompu - kridsadakorn.cha@biotec.or.th; Chumpol Ngamphiw - chumpol.nga@biotec.or.th; Uttapong Ruangrit - uttapong.rua@biotec.or.th; Juntima Chanprasert - juntima.cha@biotec.or.th; Anunchai Assawamakin - assawamakin.a@gmail.com; Sissades Tongsima* - sissades@biotec.or.th

* Corresponding author

Published: 14 August 2007

BMC Genomics 2007, 8:275 doi:10.1186/1471-2164-8-275

This article is available from: http://www.biomedcentral.com/1471-2164/8/275

© 2007 Wangkumhang et al; licensee BioMed Central Ltd.

This is an Open Access article distributed under the terms of the Creative Commons Attribution License (<u>http://creativecommons.org/licenses/by/2.0</u>), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

Received: 26 April 2007 Accepted: 14 August 2007

Abstract

Background: Allele-specific (AS) Polymerase Chain Reaction is a convenient and inexpensive method for genotyping Single Nucleotide Polymorphisms (SNPs) and mutations. It is applied in many recent studies including population genetics, molecular genetics and pharmacogenomics. Using known AS primer design tools to create primers leads to cumbersome process to inexperience users since information about SNP/mutation must be acquired from public databases prior to the design. Furthermore, most of these tools do not offer the mismatch enhancement to designed primers. The available web applications do not provide user-friendly graphical input interface and intuitive visualization of their primer results.

Results: This work presents a web-based AS primer design application called WASP. This tool can efficiently design AS primers for human SNPs as well as mutations. To assist scientists with collecting necessary information about target polymorphisms, this tool provides a local SNP database containing over 10 million SNPs of various populations from public domain databases, namely NCBI dbSNP, HapMap and JSNP respectively. This database is tightly integrated with the tool so that users can perform the design for existing SNPs without going off the site. To guarantee specificity of AS primers, the proposed system incorporates a primer specificity enhancement technique widely used in experiment protocol. In particular, WASP makes use of different destabilizing effects by introducing one deliberate 'mismatch' at the penultimate (second to last of the 3'-end) base of AS primers to improve the resulting AS primers. Furthermore, WASP offers graphical user interface through *scalable vector graphic* (SVG) draw that allow users to select SNPs and graphically visualize designed primers and their conditions.

Conclusion: WASP offers a tool for designing AS primers for both SNPs and mutations. By integrating the database for known SNPs (using gene ID or *rs* number), this tool facilitates the awkward process of getting flanking sequences and other related information from public SNP databases. It takes into account the underlying destabilizing effect to ensure the effectiveness of designed primers. With user-friendly SVG interface, WASP intuitively presents resulting designed primers, which assist users to export or to make further adjustment to the design. This software can be freely accessed at <u>http://bioinfo.biotec.or.th/WASP</u>.

Background

One of the results from human genome project is the use of human variation data called single nucleotide polymorphism (SNP). Such polymorphisms including rare ones such as mutations reflect how we have different genetic responses to the environment as well as predispositions to certain diseases and varying effects to drug treatments. These variations, therefore, play a crucial part in many studies including population genetics, molecular genetics and pharmacogenomics. Consequently, to capture the underlying variations, many genotyping technologies have been proposed, such as restriction enzyme assay, hybridization assay, ligation assay, invasive cleavage assay and allele-specific polymerase chain reaction assay. Each of these techniques has both advantages and disadvantages. In this work, we focus on a simple and cost-effective protocol, namely the allele-specific polymerase chain reaction or AS-PCR genotyping technique. This assay protocol were adopted by many recent studies, for instance, the detection of common alleles of cytochrome P450 that mediates different level of drug metabolisms [1], the screening of mutations responsible for beta-thalassemia [2], and the evaluation of using a set of SNPs as genetic markers of a disease [3].

In literature, AS-PCR is also known as amplification refractory mutation system (ARMS) [4]. This technique is a quick and dependable genotyping protocol that requires minimal instruments found in most laboratories. It is based on the extension of primer only when its 3'-end is a perfect complement to the allele present in the input sample. Thus, if a single base polymorphism occurs, the genotyping results can be observed by simply comparing the length of PCR products. Several automated bioinformatics tools are available for specifically designing AS primer, such as, Tetra-primer [5], Primo SNP [6], and a commercial software Visual-OMP [7]. However, none of the available tools offer an integrated input system that users can retrieve markers to be genotyped from within the program framework. Furthermore, the resulting AS primers from these programs may not be specific enough; that is only one mismatch at the terminal base may not discriminate the amplification results from wild and mutant types since the extension of the allele-nonspecific primer may continues if terminal mismatching has weak-destabilizing effect [8].

To address these difficulties a SNP database and intuitive graphical interface should be integrated with the AS primer design tool while primer-destabilizing condition must be considered to ensure effectiveness of primer results. Consequently, a novel tool, WASP (Web-based Allele-Specific PCR Primer designing tool), is presented in this paper. This tool automatically generates well-calibrated discrimination conditions for AS-PCR assay. By adopting the protocol, our proposed tool can increase the AS primer specificity during PCR reaction by creating a mismatch on the *penultimate* (second to the terminal) base of the AS primer. WASP utilizes a local SNP database, comprising SNPs from dbSNP [9], HapMap [10], and JSNP [11], to search and retrieve input SNP information from within the tool. Therefore, this eliminates extra steps of gathering SNPs from different websites and copying/ pasting to a primer-designing tool. Through a unique SVG graphical user interface, WASP offers a primer design for multiple SNPs that store in public databases. Note also that WASP calculates AS-PCR assay for genotyping SNPs, which are not in this database by providing flanking sequences as well as the target SNP(s) to the program. Optionally to optimize PCR experiment, standard PCR conditions can be modified through the primer-parameters. If design conditions are met, AS primers and their common oligo-primers will be graphically visualized or otherwise displayed the reasons why they cannot be designed. Once resulting AS primer set is obtained, these primers can also be verified for their global uniqueness by running through an in silico PCR [12].

Implementation

WASP is a web application constructed using Ruby on Rails Framework. The application is integrated into the local SNP database (shown in Figure 1) running MySQL database server. This database frequently collects public SNP information and related reference data from public SNP databases. At the time of writing current version of these databases are dbSNP build 126, HapMap public release 20 and JSNP release 28. The computing core runs within UNIX environment comprising four modules: 1) the input module, 2) the parsing module, 3) the AS-PCR assay analysis module and 4) the graphic display module.

I) Input module

The purpose of this module is to read in target SNP(s) or mutation(s) that users want to construct an AS-PCR assay. WASP uses two alternative input formats, namely keyword and user-defined input. First, the keyword search is designed for looking up related SNP information such as possible alleles and their corresponding flanking sequences. The allowed keywords can be Gene name and reference SNP rs number. The input module creates a set of SQL queries to retrieve underlying information such as the 5' and 3' flanking sequences and the target alleles of SNP; then pass the resulting information to the parsing module to further prepare inputs for Primer3 program (see AS-PCR assay analysis module). For visualization and selection of SNPs to be genotyped, an SVG-based graphical view of SNP locations will also be displayed once a gene is selected.



Figure I

WASP computational flow diagram. (1) the input module, (2) the parsing module, (3) the AS-PCR assay analysis module and (4) the graphic display module.

The second input format is constructed to facilitate the AS primer design for novel SNPs and/or mutations. The target variations and their corresponding 5' and 3' flanking sequences must be entered to the given text box or upload a file containing the information to the server. Instead of constructing SQL queries described previously, the input data will be converted to the internal forms by the parsing module and later sent to the AS primer analysis module. For example, Figure 2 demonstrates two possible sample inputs that WASP would except to calculate the corresponding AS primer(s). The input format for this section comprises 5' and 3' flanking sequences or a target variation, which is described in IUPAC (Figure 2a) or bracketing format (Figure 2b). Note that for both inputs (querying from database and manually entering), their primer conditions can be modified such as primer length, melting temperature (Tm), GC content for primer oligos and other standard primer design parameters.

2) Parsing module

This module is responsible for transforming the underlying variation information and their primer conditions to PCR template sequences and primer selection criteria that are used in the Primer3 program [13]. Both *wild* and mutant type templates are generated separately, however, they must share the same primer selection criteria in order to be performed in single PCR reaction. The main contribution of this module is the preparation of an additional deliberate mismatch at the penultimate base of AS primers to increase the reaction specificity [8]. For instance, if there exists a weak destabilizing at the terminal mismatch, it is likely that such a mismatch is not enough to stop the extension by Taq polymerase. To increase the discriminative power, hence, an additional mismatch with a strong destabilizing effect should be introduced at the penultimate base of the AS primers. Figure 3 demonstrates the decision criteria, based on the destabilizing conditions, in which the underlying additional mismatch should be added to AS primer. The areas shaded in red and yellow respectively represent the strong and weak destabilizing effect of any two nucleotides while blue ones imply the medium effect. Starting from the center of this figure (3' end of AS primer) and working outward, we can decide which base (the outermost base of this diagram) to be placed on the penultimate position of the designed AS primer. An example situation showed under the diagram, AS primer was designed for genotyping C>T SNP at +2850 locus of the CYP2D6 gene [14]. The polymorphisms, C

>SNP_1 ...CCCATCCTTTTCGAGTTGAAGTT ATTTTGAGAACACATCCCCTAAGW GTTTATCTTGGTTATTTGC...

(a)

>SNP_1 ...CCCATCCTTTTCGAGTTGAAGTT ATTTTGAGAACACATCCCCTAAG[A/T] GTTTATCTTGGTTATTTGC...

(b)

Figure 2 WASP sample input. (a) IUPAC ambiguity code and (b) Bracketing code format

(wild-type) and T (mutant), are displayed in the box on the template sequence (3' to 5' direction). On top of each template, primer extension direction is indicated by \rightarrow whereas × marks a non-extended scenario. Note that nucleotides in the AS primers that are not complementary to the target appear afloat from the primer oligos. Basically, four binding scenarios could happen in a single AS-PCR experiment: WP-WT, WP-MT, MP-WT and MP-MT. To design the AS primers the above polymorphism, wild type and mutant AS primers must have respectively bases G and A at the 3' terminus; thus gel electrophoresis bands can be observed when successful hybridization occurs.

In order to enhance reaction specificity, a deliberate mismatch is introduced at the penultimate base. Using the above diagram, we first consider wild type primer with the case WP-MT in which T/G (^①) presents a weak-destabilizing mismatch (yellow). We then look at the penultimate location to observe the base on the template sequence, which is G (^②) in this case. According to the diagram only nucleotide A (^③) should be put next to the terminal base in both AS primers to ensure the stronger destabilizing (red). Note here that the same additional mismatch is also introduced to mutant primer. If the terminal base is not mismatching, the additional mismatch on the penultimate base (^④ and ^⑤) of the primer will not prevent the extension.

3) AS-PCR assay analysis module

Three oligo primers are to be produced by the parsing module, namely wild-type, mutant and common primers. The AS-PCR assay will be conducted in two parallel experiments: one is "wild+common" primer experiment (WC) and the other is "mutant+common" primer experiment (MC). The results will be interpreted from each experimental result. Three scenarios are to be detected, which are homozygous wild type, heterozygous and homozygous variant type respectively. A band from WC but none from MC tube entails homozygous wild type. On the other hand, single band on MC but none from WC tube implies homozygous variant type. If this kind of band appears on both WC and MC tubes, being heterozygous will be inferred.

With the above AS-PCR assay, the core engine of this module is the Primer3 program with some modification to make Primer3 fit the needs of AS primer designing. Two criteria must be considered for this AS-PCR assay. First, AS primers must be used with their common primer in the same reaction. Secondly, the reaction conditions should promote the discrimination of wild and mutant reactions. In this module, optimal primer candidates will be computed. Both wild and mutant primers are examined if they follow given AS-PCR assay criteria. Once satisfied, each primer pair will be aligned against the human genome reference sequence stored in the local database using in silico PCR. This technique utilizes BLAT (BLAST-Like Alignment Tool) to rapidly search for all possible PCR products in the target genome that complement with both directions of PCR primers. In our case, the input of in silico PCR is a pair of optimal primer candidates and the output returns PCR products if any. This feature ensures that the global uniqueness of a primer pair can be achieved.

4) Graphic display module

The graphic display module is responsible for graphically displaying the query SNPs as well as the designed primer(s). We adopted the scalable vector graphic (SVG) display [15] that allows users to interact with high quality vector image that can be scaled and move around without viewing degradation due to losses in image compression. The construction of such images is performed "on-demand", i.e., no images are pre-calculated and stored. First, for SNP location map (see Figure 4c), relative locations on each gene are computed upon the query to the local SNP database. On SVG enabled web browsers such as Firefox and Opera SVG images can be displayed without extra installation while all other browsers SVG plug-in module must be installed.

Second, for displaying AS primers results, the display module waits for results from AS-PCR assay analysis module to be completed. Then the primer information such as



Figure 3 Graphical view of destabilizing effect adapted from ARMS done by Little [7].

orientation, location, and other calculated conditions are converted to SVG image. Here, target SNPs/mutations and the resulting AS primers will be visualized in graphical objects. When moving a mouse over each object, there will be a popup window whose contents show related primer conditions. If no AS primer is reported, only input sequence and original primer parameter conditions will be presented. In this case, some primer selection criteria should be optimized such as varying primer length, Tm, GC content as well as increasing the 5'and 3' flanking sequence length. Figure 4 demonstrates SVG graphical user interfaces described previously using AS primer design for SNPs in CYP2D6. This example demonstrates the case when multiple SNPs on the gene can be selected for AS primer design (Figure 4e). Figure 5 presents primer report (Figure 5a–5b) as well as the SVG output representation (Figure 5c). The results are collectively reported in one continuous page if there are multiple SNPs to be designed. Such results can directly be saved as a text file for future experiments.

	ſ	1							_				1					
		۷	Veb-bo	ASP ased Allel	le Sp	ecific Prime	lesigning to I			Home	Abo	ut Us Contact Help					-	
					Sect	ion A: Allele-s	/	Section	n A: Alle	le-spec	ific	primers of	SNPs in	the E	Databa	Ise		
					Searc	h by Genero		Search b	Gene	Name 🛟	for	CYP2D6		in 💽	dbSNP	•		
					Contin	ue Section C if you wish to	adjust primer se	lection parameter	ers then click subm	it below							-	
		Ass File Mar	tion B: Us ay Name (O; Uploading: hual Input:	ser-defined Inp ptional): Choose File no I	file selecto	d		Section C: I Primer Lengtl GC% Optimal: Tm Optimal: [Max Tm Diff: Max Self Com Max 3' Self Co Max polyX in Number prim Primer Filterii	Primer3 Parai h Optimal: 20 50 Min: 4 100 plementarity: complementarity primer: 3 er pairs to retu	Min: 18 20 Max: 81 0 Max: 65 8	Max 5	t c 36					a	
1. Gene: CYP	200	5							ng.	•				1				
Difficial symbol <u>YP2D6</u> Full name : cy Refseq status : REVIEWED Polymorphism Insert on : Co	ytoc	G hrom S	e P450, fa Disorder : ensitivity	1565 mily 2, subfami Parkinsonis	ily D, po m, susc	lypeptide 6 eptibility to (1) De	brisoquine	OMIM Id : Alias : Cf P450-DB1, Chromoso	: <u>MIM : 124</u> PD6, CYP2D, 0 , P450C2D ome : 22	030 G CYP2D@, CYP C	ene 2DL1 ytog	type: protein-cod , MGC120389, MGC enetic band: 22c	ing 120390, 113.1				b	
sSNP : 136 (Va	alida	ted 2	8 SNPs)	Hapmap : 2	JSNP :	<u>3</u>	•											
fSNP on Gene I	ID :	1565																
3' Reverse								Wholem	nap (4383 bp)						5'		
21912600		2	1913400	2191	14200	21915000	0	21915800)	21916600		21917400	219182	00	2191	NT_011520.	11	
40852031 RefSNP		4	0852831	4085	53631	40854431		40855231		40856031		40856831	408576.	51	4085	6431 Chr 23	2	C
	i		-				II dente i i				-		111					C
				locus i	region	intron	exon	1	new	nonsync	on	validated	verified					
					c	Gene CYP2D6 on mRI	NA NM_0001	06.4 conti	g NT_011520.	11 Chr 22 (h	ave	9 exons)	Сору	right 2007	BIOTEC. A	II Rights Reser	ved.	
SVG help	Γ.																_	
		Ref	erence SNP	(dbSNP) Repor	nt								Bael					
		Mole	cule Type		genomi	:	db	SNP ID : rs161	947									
	1	Varia	tion Class		snp A/G												-	
		Valid	ation Status	s (<u>Helo</u>)	N												d	
a d	7	2016	047 LalieleR	te 121 totali ene 1	241		Fla	inking Sequer	nces								u	
A C	C A A	CTGT	TTCCCA GAT	GGGCTCA CGCTGC	ACAT CC	GGATGTAG GATCATGAG	C AGGAGGCC	20										
~	C A	A/G CAGG	ттетса тел	TTGAAGC TGCTCT	CAGG GT	тессетта асстанася	G GGCCGAGAG	ac .										
A		ATAC	COUVA CAG	ANCOUNT TAUCCO	SUAA AT	ATTETOCAC	TOTCAGCCO								_	i .		
2	A1 S1	Avera	Selected	Chin 10		Volidation (Male)	Senomic Da	Minor	Position	Chr	-	Tr Des Tress	anscription Di	Codon	AA			
		Stand	đ	SNP ID	туре	validation (<u>Heip</u>)	ouserved	Allele	on Contig	Position	str	Pxn Type	Aut change	Position	Position			
		~	 ✓ 	rs1058172	snp	F	C/T	C: 0.004	21913996	40853427	F	coding-nonsynon	H/R	2	365			
		22	2	rs3915951	snp	F	G/T	G: 0.380	21914149	40853580	R	coding-nonsynon	L/R	2	329			
		22	2	rs1058170	snp	F	C/G	C:0.065	21914443	40853874	R	coding-nonsynon	G/A	2	300			
	H		2	rs16947	snp	N	A/G		21914456	40853887	F	coding-nonsynon	R/C	1	296		0	
		Chr	2	rs17002852	snp	CF	C/T	C:0.440	21914836	40854267	R	coding-synon	H/H	3	232		e	
		22	2	rs17002853	snp	F	C/T	C:0.340	21914840	40854271	R	coding-nonsynon	P/L	2	231		-	
4		22	2	rs4996602	snp	N	A/G	A:0.003	21915058	40854489	F	intron						
	٦	22		rs2267447	snp	CF2	C/T	C:0.413	21915209	40854640	F	intron	610		174			
				rs11569739	snp	F 2	A/G	A : 0.419	21915437	40854868	F	coding-synon	G/G	3	1/6			
			-	1011000728	snp	F	70	A 10.169	21913459	40034090	R	coung-nonsynon	E/G	2	109			
				rs1058164	Son.	F	CIG	C:0340	21915645	40855076	P	coding-suppor	VA	3	136			

Figure 4

WASP SVG graphical user input interface. (a) WASP first page containing three input sections: 1) keyword input for local SNP db, 2) user-defined input and 3) AS primer parameter adjustment. (b) CYP2D6 gene output resulting from using CYP2D6 as a keyword search. The general information including gene locus, gene type, OMIM Id and related disease are reported here. Polymorphism related information on CYP2D6 gene is also shown e.g. number of SNPs in CYP2D6 gene according to dbSNP, HapMap and JSNP. (c) Gene structure view of CYP2D6 displayed as SVG image. The top ruler displays positions and a direction of the gene relative to the reference Contig and chromosome respectively. (d) Detailed information of each SNP including flanking DNA sequences, observed alleles, allele frequency, position on Contig/chromosome, amino acid change and its position relative to the transcription. (e) Selectable SNP list for AS primer design. Using the previously described information, users can select a set of SNP from this table for designing AS primers.

leverse						Wholemap	(4383 bp)							5'	
)	21913400	2	1914200	21915000		21915800	2	19166	500 219	17400	219182	00	21919	NT_011520.	11
1	40852831	4	0853631	40854431		40855231	4(08560	031 408	56831	408576	31	40858	431	
														Chr 2	2
		l ilini			llinen 14			_							
	-	- locu	is region – in	ntron	exon	nev	v	no	nsynon valio	lated	verified				
			Gene CYP	2D6 on mRN	A NM_00010	6.4 contig N	T_011520.11	Chr	22 (have 9 exons)						
	a		0110 M D				111 - 11 - 11 - 11 - 11 - 11 - 11 - 11				Copy	yright 2007	BIOTEC. All	Rights Reser	rved.
help Zoom ii	n : Ctrl + Mouse, Zoom	out : Ctrl	+ Shift + Mouse, Pan	window : Alt	+ Mouse, Sho	v menu : Right	click (On Macint	losh re	eplace Ctrl with Apple-	Key)					
				Genomic Dat	ta				Tra	inscription Da	ata		1		
: Res	ult	7.000	Validation (Male)	Observed	Minor	Position	Chr	6 4-	Die Dies		Codon	AA			
	SNP ID	туре	validation (Help)	Observed	Allele	on Contig	Position	Str	Exn Type	AA Change	Position	Position			
	rs2982054	snp	N	С/Т	C:0.004	21913996	40853427	F	coding-nonsynon	H/R	2	380			
**	rs1058172	snp	F	A/G	A:0.280	21914041	40853472	R	coding-nonsynon	H/R	2	365			
	rs3915951	snp	F	G/T	G:0.380	21914149	40853580	R	coding-nonsynon	L/R	2	329			
\vee	rs1058170	snp	F	C/G	C:0.065	21914443	40853874	R	coding-nonsynon	G/A	2	300		h	
	rs16947	snp	Ν	A/G		21914456	40853887	F	coding-nonsynon	R/C	1	296		N	
*	rs17002852	snp	CF	С/Т	C:0.440	21914836	40854267	R	coding-synon	н/н	3	232			
*	rs17002853	snp	F	С/Т	C:0.340	21914840	40854271	R	coding-nonsynon	P/L	2	231			
Ň	rs4996602	snp	N	A/G	A:0.003	21915058	40854489	F	intron						
×	rs2267447	snp	CF2	с/т	C:0.413	21915209	40854640	F	intron						
*	rs2267448	snp	F 2	A/G	A: 0.419	21915437	40854868	F	coding-synon	G/G	3	176			
			Opt:	:20 Min:18 Ma	exitectacies	CAGGTGGGCCCAL	TETAGGAACCETGG	CCAC	CTAGTCCTCAATG			ATCOARGIACI	т		
	GC% TM		Opt: Opt: Opt:	20 Min: 18 Ma :50.0 Min: 20. :55.0 Min: 40.	0 Max:85.0 0 Max:65.0	CAGGTGGGCCCAC	TCTAGGAACCCTGG	ICCACI	CEAGECCECAAEG		8.0	AICOARGIACI	T -		
	GC% TM Max Tm Diff Max 3' Complement	ntarity	Opt: Opt: Opt: 100. 3.0	20 Min:18 Ma 50.0 Min:20. 55.0 Min:40.	ax:36 0 Max:85.0 0 Max:65.0	CAGGTGGGCCCAC	TCTAGGAACCCTGG	Max S	Self Complementarity		8.0 3	ATCOMOTACI	T		
	GC% TM Max Tm Diff Max 3' Complement Number Primer Ro	ntarity eturn	Opt: Opt: Opt: 100. 3.0 S	:20 Min:18 Ma :50.0 Min:20.1 :55.0 Min:40.1	CATGUTACTU ax:36 0 Max:85.0 0 Max:65.0	CAUGUCCAU	Mismatch	Max S Max S	Self Complementarity ax PolyX in Primer he penultimate primer	position	8.0 3 Yes	ALCOARTAC	T		
	GC% TM Max Tm Diff Max 3' Complemen Number Primer R In-silico PCR Filte	ntarity eturn ering	Opt: Opt: Opt: 100. 3.0 S No	20 Min: 18 Ma 50.0 Min: 20.0 55.0 Min: 40.0	CATGOCTACTUC ax:36 0 Max:85.0 0 Max:65.0	CAUGUCCAU	Mismatch	Max S Mat I Mi at th	Self Complementarity ax PolyX in Primer he penultimate primer lispriming Filtering	position	8.0 3 Yes No				
rs3915951	GC% TM Max Tm Diff Max 3' Complement Number Primer R In-silico PCR Filte	ntarity eturn ering	Opt: Opt: 000 3.0 5 No	20 Min:18 Ma 50.0 Min:20. 55.0 Min:40.	erroutroto ax:36 0 Max:85.0 0 Max:65.0		Mismatch	Max S Max S Ma at th	Self Complementarity ax PolyX in Primer he penultimate primer Ispriming Filtering	position	8.0 3 Yes No				
rs3915951 91	GC% TM Max Tm Diff Max 3' Complemen Number Primer Ri In-silico PCR Filto 121 15	ntarity eturn ering	0pt: 0pt: 100. 3.0 5 No 181 211	20 Min: 18 Ma 50.0 Min: 20.0 55.0 Min: 40.0 0	271	301	Mismatch 	Max S Ma at th M	Self Complementarity as PolyX in Primer he penultimate primer lispriming Filtering 391 421	position 451	8.0 3 Yes No 481	511	T -		
rs3915951 91	GC% TM Max Tm Diff Max 3' Complement Number Primer R In-silico PCR Filtz 121 15	ntarity eturn ering il	0pt: 0pt: 0pt: 100. 3.0 5 No 181 211	20 Min: 18 Mi 50.0 Min: 20.0 55.0 Min: 40.0 0 241	271	301 7s391595	Mismatch	Max S Max S Mat th M	Self Complementarity ax PolyX in Primer he penultimate primer lispriming Filtering 391 421	position 451	8.0 3 Yes No	511		1	
- rs3915951 91 I	GC% TM Max 3' Complement Number Primer R In-silico PCR Filta 121 15	ntarity eturn ering	Opt: Opt: 100. 3.0 5 No 181 211	20 Min: 18 Mi 50.0 Min: 20.0 55.0 Min: 40.0 0	271	30 rs39156 Alte: Cl Position: 3	Mismatch	Max S Ma at th	Self Complementarity ax PolyX in Primer he penultimate primer lispriming Filtering	position 451	8.0 3 Yes No 481	511			
91 91 0ligo1 [Prodi	GC% TM Max Tm Diff Max 3° Complement Number Primer R In-silico PCR Filtx 121 15 1	ntarity eturn ering	0 0pt: 0 0pt: 100, 3,0 8 181 211 181 211	20 Min: 20 Min: 20. : 50.0 Min: 20. : 55.0 Min: 40. 0 241	281 281	304 ra391569 Poaltion: 3	Mismatch 	Max S Ma at th M	Self Complementarity ax PolyX in Primer he penultimate primer Ispriming Filtering	position 451	8.0 3 Yes No 481	511 			
rs3915951 91 1 Oligo1 [Produ	GC% TM Max Tr Diff Max 3' Complemen Number Primer R In-silico PCR Flitt 121 15 121 15 121 15	ntarity eturn ering il Wie do Tm	00t: 00t: 100. 3.0 8 No 181 211 191 211 40t forest from t Apple forest from t debatter f debatter f debatter f debatter f debatter f debatter	20 Min: 18 Min: 20. : 50.0 Min: 20. : 55.0 Min: 40. : 0 241 : 55.0 Min: 40. : 0 : 55.0 Min: 40. : 55.0 Min: 40. : 55.0 Min: 40. : 55.0 Min: 40. : 50.0 Min: 20. : 50.0 Min: 40. : 50	271 281 281	30 Adject of Position 3	Mismatch 	Max S Ma at th M	Self Complementarity ax PolyX in Primer he penultimate primer ispriming Filtering	position 451 460	8.0 3 Yes No 481	511 			
rs3915951 91 1 Oligo1 [Prodi	GC% TM Max Tm Diff Max 32 Compleme Number Primer R In-silico PCR Filt 121 15 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	ntarity eturn ering ii Wite Goi Tim Mite Goi	Opt: Opt: Opt: Opt: 100. 3.0 No 181 211 Ryst Forward Primer 9* No 184 211 Ryst Forward Primer 9* Rest Optimer 9* Rest Optimer 9* eff Optimer Primer 9* Rest Optimer 9* eff Optimer Primer 9* eff Optimer 9* eff Optimer 9*	20 Min 13 B 420 Min 13 Min 140	271 271 281 281 281	30 Adecid	Mismatch 	Max S Max Max Max Max	Self Complementarity ax PolyX in Primer he penultimate primer lispriming Filtering	position 451 460	8.0 3 Yes No 481 Common Reve CATGTCACCC Tmis4.90, GC	511 	10 sp)	2	(
	GC% TM Max Tm Diff Max 3 ⁻ Complement Number Primer R In-silico PCR Filts 121 15 1 121 15 1 100 PCR Filts 100 PCR	ntarity eturn ering ii Wite GG Tm	Opt: Opt: Opt: Opt: 100. 3.0 No 181 211 hps Forward Primer 9' 0.00, CO-40.0, Set Anyot 43.07, CO-40.0, Set Anyot	241 241 50.0 Min: 20.0 Min: 20.0 55.0 Min: 40.0 0 241 1 10 bp) 6.0, Set End-0.00 0 bp) 0 bp) 0 bp) 0 bp) 0 bp)	271 281 281 281	30 rs3158 Atec C	Mismatch 	Max 1 Max 1 Ma at th M	Self Complementarity ax PolyX in Primer he penultimate primer lispriming Filtering	position 451 460	8.0 3 Yes No 481	511 511 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	P P P P P P P P P P P P P P P P P P P	1	
S3915951 91 	GC% TM Max Tm Diff Max 3 Complemen Number Print In-silico PCR Filto 121 15 121 15 15 15 15 15 15 15 15 15 15 15 15 15 1	ntarity eturn ering ii ii Wite Goi Tim Mut Goi Tim Wite Goi Tim Wite Goi Tim	Opt:	241 241 241 260 or 100 or 1	271 281 281 281 281	30 rs31562 Atel of Position 2	Mismatch 	Max S Ma at th M	Self Complementarity ax PolyX in Primer he penultimate primer ispriming Filtering	451	8.0 3 Yes No 481	511 TRE Prime S DOATGTOATA (TRE Prime S) DOATGTOATA (TRE Prime S)	10 bp)		
S3915951 91 	GC% Comparison of the second s	ntarity eturn ering il Wid Gd Tm Wid Gd Tm	Opt: Opt: Opt: Opt: Opt: 100. 3.0 No 181 211 + 40.0 24.0 24.0 24.0 24.0 24.0 24.0 25.0 21.0	201 Min: 38 Min: 20 Min: 38 Min: 20 Min: 38 Min: 20 Min: 38 Min: 20 Mi	271 281 281 281 281	30 Alector Carl	Mismatch 	Max 1 Max 1 Ma at th M	Self Complementarity ax PolyX in Primer he penultimate primer ispriming Filtering	position 451 460	8.0 3 Yes No 481 Common Reve ArrorCACCC Tm-64.93, GC	511 	10 top) 10		
Oligo1 (Prodi	GC% TM Max Tm Diff Max 3 ² Complement Number Primers In-silico PCR Filt 121 15 121 15 121 15 121 15 121 15 121 15 121 15 121 15	itarity eturn ering ii Multi GG Tm Tm Will GG GG	Opt: Opt: Opt: Opt: Opt: Opt: Itel:	20 Min: 38 Min: 20. 50.0 Min: 12 Min: 20. 55.0 Min: 20. 55.0 Min: 20. 50.0 Min: 20. 241 1 0 Min: 20. 50.0 Min: 20. 50	271 281 281 281 282 282	30 Address of the second	Mismatch 	Max 1 Max 1 Ma	Self Complementarity ax PolyX in Primer he penultimate primer ispriming Filtering	position 451 460 480	8.0 3 Yes No 481 Common Reve Arroro.coco Tm-54.93, 00 Tm-54.93, 00	511 Tree Primer S ¹ Source Carlo Source Carlo Source Carlo Source S ¹ Source S ¹ Sourc	9 bp) 19 bp)		
Dilgo1 (Prodi	GC% TM Max Tm Diff Max 3 ² Complement Number Prime R In-silico PCR Filt 121 15 121 15 15 15 15 15 15 15 15 15 15 15 15 15 1	ntarity eturn ering ii Mada Tim Wile Ga Tim Tim Tim	Opt:	2241 2241 235.0 Min 13 Min 140. 0 241 10 lps) 241 10 lps) 241 10 lps) 0 l	271 281 281 281 281 282 282 282	30 Addition	Mismatch	Max S Ma a at th M	Self Complementarity ax PolyX in Primer he penultimate primer Ispirming Filtering	position	8.0 3 Yes No 481 1 Common Reve GATGTCACO Timode 80, GC	511 srse Primer S [*] GOATGTCATA (Simon 27, See Arriver arse Primer S [*] GOATGTCATA (Simon 27, See Arriver Simon 27, See Arriver	19 bp) 19 color. Set End-3.	×	
Cligo1 [Prodi	GC% TM Max Tm Diff Max 3 ⁻ Compleme Number Primer R 1n-silico PCR Filt 121 15 121 15 15 15 15 15 15 15 15 15 15 15 15 15 1	ntarity eturn iii Mada Tim Mada Tim Tim Tim Tim Tim Tim	Opt: Opt: Opt: Opt: Opt: Itop:	2010 Min 28 Min 200 50.0 Min	271 271 281 281 282 282 282 282 282	30 Pastici	Mismatch 	Max 1 Max 1 Max 1 Max 1 Max 1 Max 1 Max 1	Self Complementarity ax PolyX in Primer he penultimate primer ispriming Filtering	position 451 1 460 460	8.0 3 Yes No 481 Common Reve GATGTCACCC Time54.93, GC	511 	19 hp) 19 hp) 19 hp) 19 hp) 19 hp) 19 hp)		
Oligo1 (Prodi	GC% Compension of the second s	ntarity eturn ii Waa Goo Tmi Moo Goo Tmi Moo Goo Tmi Moo Goo Tmi Moo Goo Tmi Moo Goo Tmi Moo Goo Tmi Moo Goo Tmi Moo Goo Tmi Moo Coo Coo Tmi Moo Coo Coo Coo Coo Coo Coo Coo Coo Coo	Opt: Opt: Opt: Opt: Opt: 100. 100. 100. 100. 100. 100. 100. 100. 101. 102. 103. 104. 104. 105. 106.	241 241 1 0 0 0 0 0 0 0 0 0 0 0 0 0	271 271 281 281 282 282 282 282	30 Address Add	Mismatch 	Max S Max Ma	Self Complementarity ax Polyk in Primer he penultimate primer ispriming Filtering	position	8.0 3 Yes No 481 Common Reve Garagroad Transfer Garagroad Garagroa	511 	9 kg) 9 kg 9		
- rs3915951 91 Oligo1 (Prodi Oligo2 (Prodi	GC% TM Max Tm Diff Max 3 ⁻ Complement Number Print In-silico PCR Filts 121 15 121 15 15 15 15 15 15 15 15 15 15 15 15 15 1	ntarity eturn ii Weider Good The Model Good The The Model Good The The Model Good The Model Good The Model Good The Model Cool Cool Cool Cool Cool Cool Cool Co	Opt: Opt: Opt: Opt: Opt: 100. 3.0 No 181 211 hps:Forward Nimer F 0.00, CO-40.00, SM Any- 0.	241 241 241 241 241 25.0 Min 12 Min 20. 25.0 Min 12 Min 20. 25.0 Min 2	271 271 281 281 281 282 282 282 283	30 Another and a second s	Mismatch 	Max S	Self Complementarity ax PolyX in Primer he penultimate primer ispriming Filtering	position	8.0 3 Yes No 481 Common Reverse Anternexcel Common Reverse Common	511 	9 bg) 19 bg 19 bg		
-s3915951 91 Oligo1 (Prodi Oligo2 (Prodi	GC% Comparison of Comparison o	ntarity eturn i1 i1 i2 i1 i1 i1 i1 i1 i1 i1 i1 i1 i1 i1 i1 i1	Opt: Opt: Opt: Opt: Opt: ID:	20 Min: 12 Min: 20 Min	271 271 281 281 282 282 282 283 283	30 Able Call	Mismatch 	Max 5 Ma th M	Self Complementarity ax PolyX in Primer he penultimate primer ispriming Filtering	position 451 460 460 460 460 460	8.0 3 Ves No 481 481 Common Reserve Autor Charlos Common Reserve Common Reserve Reserve Common Reserve Reserve Common Reserve Common Reserve Reserve Common Reserve Common Reserve Reserve Common Reserve Common Reserve	511 The Prime S Source Carl (Source Carl (Sou	19 bp) 19 bp)		
S3915951 91 Oligo1 [Prodi Oligo2 [Prodi Oligo3 [Prodi	GC% UNDER CONTRACTOR OF CONTRA	starity etum aring ii Web GG GT The The Man Man The The Man Man The Man Man The Man Man The Man Man Man Man Man Man Man Man Man Man	Opt: Opt: Opt: Opt: Itio:	20 Min: 18 Min: 20. 50.0 Min: 18 Min: 20. 55.0 Min: 20. 55.0 Min: 20. 55.0 Min: 20. 50.0 Min: 20. 241 10 Min: 20. 50.0 Min: 20. 50.	271 271 281 281 282 282 282 283 283 283		Mismatch	Max S Ma at th M	Self Complementarity ax PolyX in Primer he penultimate primer ispriming Filtering	position 451 460 460 460	8.0 3 Ves No 481 4 4 Common Reveal Antorecological Common Reveal Antorecological Tim-set 49, 00 Common Reveal Tim-set 49, 00	511 THE Prime 9 THE PRIME 9	9 bp) 19 bp) 10 bp)		
S3915951 91 Oligo1 [Prodi Oligo2 [Prodi Oligo3 [Prodi	GC% TM TM Max Tm Diff Max 7 Compleme Max 7 Compleme In-selfco PCR Filt 121 15 15 15 15 15 15 15 15 15 15 15 15 15 1	starity etum ering ii Weid did Tm Tm Weid did did Tm Tm Weid did did tm Tm Midd did did tm Tm Midd did did tm Tm Midd did did tm Midd did did tm Midd did did tm Midd did did tm Midd did did did tm Midd did did did did did did did did di	Opt: Opt: Opt: Opt: ID0.	20 Min: 28 Min: 20. 550.0 Min: 20. 550.0 Min: 20. 550.0 Min: 20. 50.0 Min: 20. 50.0 Min: 20. 241. 1 1 1 1 1 1 1 1 1 1 1 1 1	271 271 281 281 282 282 282 283 283	300 Register Position 300 200 200 200 200 200 200 200	Mismatch 	Max 1 Max 1	Self Complementarity ax PolyX in Primer he penultimate primer ispriming Filtering 301 421 3 4 4	position	8.0 3 Yes No 481 481 481 481 481 481 481 481 481 481	511 1 1 1 1 1 1 1 1 1 1 1 1	19 bp) 19 bp) 19 bp) 19 bp) 19 bp) 19 bp) 19 bp) 19 bp) 19 bp) 19 bp)		
Cligo1 [Prodi Cligo2 [Prodi Cligo3 [Prodi Cligo3 [Prodi	GC% TM TM Max Tm Diff Max 2 ⁻ Complement Max 2 ⁻ Complement In-silico PCR Filte 121 15 121 15 121 15 121 15 121 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 125 125 125 125 125 125 125 125 125 1	ntarity etum i1 Mada Grim Tim Mada Grim Tim Weig Grim Tim	Opt: Opt: Opt: Opt: Opt: ID:	241 241 1 241 1 0 0 0 0 0 0 0 0 0 0 0 0 0	271 271 281 281 282 282 283 283	300 1000 1	Mismatch	Max 1 Max 1	Self Complementarity ax Polyk in Primer he penultimate primer ispriming Filtering 391 421 1	position 451 460 460 460	8.0 3 Yes No 481 4 Control block minds 80, 00 Control block minds 80, 00 Control block minds 80, 00 Control block Control block	511 	9 kg) 19 kg) 19 kg) 19 kg) 19 kg) 19 kg) 19 kg) 19 kg) 10 kg 10 kg 1		
Cligo1 (Prodi Oligo2 (Prodi Oligo3 (Prodi Oligo4 (Prodi	GC% TM GC% TM Max Tm Diff Max 3 ² Complement In-silico PCR Filte 11- 121 15 121 15 121 121 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 122 15 12 12 15 12 12 15 15 12 15 15 15 15 15 15 15 15 15 15 15 15 15	ntarity etum i1 Mada Grim Tim Mada Grim Tim Tim Tim Tim Tim Tim	Opt: Opt: Opt: Opt: Opt: Opt: Ital	241 241 241 1 260, 5ef End-200 200, 5ef End	271 271 281 281 281 282 282 282 283 283	30 733159 733159 74	Mismatch Mismatch 331 361 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	Max 1 Max 1	Self Complementanity ax PolyX in Primer he penultimate primer ispriming Filtering 391 421 4 421 4 5 6 6 6 6 6 6 6 6 7 6 7 8 7 8 7 8 7 8 7 8	position	8.0 3 Yes No 461 Common Revenue R	511 	9 by) 9 by) 9 by) 9 by) 19		
Cligo 2 (Prodi Cligo 3 (Prodi Cligo 4 (Prodi Cligo 4 (Prodi Cligo 5 (Prodi	GC% TM Max Tm Diff Max 3 ² Complement In-silico PCR Filte 121 15 15 15 15 15 15 15 15 15 15 15 15 15 1	ntarity eturn ii Www.communication final www.communication final final www.communication final f	Opt: Opt: Opt: Opt: Opt: ID:	241 241 241 1 260, Sef End-200 20,	271 271 281 281 282 282 282 283 283	30 7391682 749162 7	Mismatch Mismatch 331 361 1 1 1 1 1 1 1 1 1 1 1 1 1	Max 1 Max 1 Max 2 Max 2	Self Complementantly ax PolyX in Primer he penultimate primer ispriming Filtering 301 421 41 500 421 500 420 500 4200 500 4000 500 400000000000000000000000	position 451 450 460 460	8.0 3 Ves No 481 Control Network ANOTOCACO Control Network Control Network Con	511 The Primer 9' SGATGTCATA (The Primer 9' The Primer 9' Th	9 bp) 19 bp) 10 bp)		
Cligo 1 (Prodi Cligo 2 (Prodi Cligo 3 (Prodi Cligo 4 (Prodi Cligo 4 (Prodi Cligo 5 (Prodi	GC% TM Max Tm Diff Max 37 Complement Mamber Minner R In-silico PCR Filte 121 15 15 15 15 15 15 15 15 15 15 15 15 15 1	ntarity eturn ii Multi Multi Tim Tim Tim Tim Tim Tim Tim Tim Tim Ti	Opt: Opt: Opt: Opt: Itil:	20 Min: 18 Min: 40. 50.0 Min: 12 Min: 20. 55.0 Min: 20. 55.0 Min: 20. 55.0 Min: 40. 0	271 271 281 281 281 282 282 283 283 283	300 700 700 700 700 700 700 700	Mismatch Mis	Max 1 Max 1 Max 2 Max 2	CENSTCECEDARG Self Complementantly ax PolyX in Primer he penultimate primer ispriming Filtering 391 421 391	position	8.0 3 Ves No 481 481 481 481 481 481 481 481 481 481	511 THE Prime V SEATOTCATA (Seator State of the seator	9 bp) 9 bp)		
Sign 1 (Produ Digo 2 (Produ Digo 3 (Produ Digo 4 (Produ Digo 5 (Produ 3 C (Produ	GC% TM TM Max Tm Diff Max 27 Complement In-selfcore Primer R International Prime	ntartty eturn sring 11 Wood Green Green Tmm Mac Green Green Tmm Mac Green Mac Green Ma	Opt: Opt: Opt: Opt: ID0.	2010 Min: 28 Min: 200 550.0 Min: 200 550.0 Min: 200 550.0 Min: 200 550.0 Min: 200 241 1 1 1 1 1 1 1 1 1 1 1 1 1	271 271 281 281 282 282 282 283 283	300 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	Mismatch Mismatch Mismatch	B Primer 2	Self Complementarity ax PolyX in Primer he penultimate primer ispriming Fittering 301 421 301 421 300 421 3000000000000000000000000000000000000	position	8.0 3 Ves No 481 481 481 481 481 481 487 647 647 647 647 647 647 647 647 647 64	511 1 1 1 1 1 1 1 1 1 1 1 1			

Figure 5

Graphical output of designed AS primers. (a) SVG graphical view of selected SNPs shown by blue dots on the SNPs. (b) Table indicates if the selected SNPs have their AS primers successfully designed. (c) By clicking on the checked row in the table, the SVG view of designed AS primers will be displayed. Top 5 oligo primers are displayed for each design where green and red arrows (or) show wild and mutant AS primers respectively. The blue arrow (\leftarrow or \rightarrow) indicates the common primer. Sizes of the primers are displayed above the arrows. Upon mouse-over action on each primer, the detailed information will be presented in the box next to it.

results and Discussion

Integrating SNP database to WASP makes it possible to display SNP related information to users (as seen in Figure 4b-4d); hence we can use such information to decide which SNP in a given gene to be designed for AS primers. The proposed tool can design more effective AS primers than other available software (Tetra-primer, Primo SNP and Visual-OMP) since WASP checks the uniqueness of the primer using *in silico* PCR during the design. Moreover, by incorporating ARMS enhancement protocol to the software, such AS primer results can be even more effective as this technique has been adopted by several recent studies. Note that other existing tools only produce mere primers but not the counterpart "common primer", which needs to have very similar condition to the AS primers. WASP also offers SVG graphical user interface for both SNP selection and AS primer viewing. For SNP selection, WASP allows users to visualize SNPs from a gene structure view (showing intron, exon, 5' and 3' UTR regions of CYP2D6 gene in Figures 4c and 5a). Users need not to construct the input information by accessing other public SNP databases. Finally, unlike other tools, WASP graphically displays location of top 5 output primers as well as their common primers. This kind of outputs should help users tailor their final primers to fit the AS-PCR assay.

Conclusion

We presented WASP, a web-based allele specific primer design application. The tool conveniently assists scientists in getting SNP information from local SNP database, mirrored from all major public databases, in order to design AS primers for existing SNPs. While current AS-PCR tools such as Tetra-primer, Primo SNP and Visual-OMP, require manual entering SNP information, WASP simplifies the data preparation step by graphically displaying available SNPs to be genotyped and transparently send such information to be computed by Primer3 program. This tool also offers a conventional input for either confirming the information in the database or designing AS primers for novel SNP/mutation.

To improve the specificity of PCR reactions, WASP deploys an *in silico* PCR by searching if a given primer is unique in the genome. Furthermore, WASP offers a missing feature in other primer design tool to improve the effectiveness of the resulting AS primers; that is the introduction of a proper mismatch to the penultimate nucleotide next to the allele specific site in the primer. SVG visualization is integrated into WASP so that primer results can be displayed interactively and intuitively; thus helping users to adjust primer parameters for the negative AS primer outputs. For positive AS primer results, users can export the primer conditions as a text file.

Availability and Requirements

WASP sever is available for public access at <u>http://bio</u> info.biotec.or.th/WASP

WASP requires SVG enabled web browser Opera and Firefox, which are freely available on almost all major Operating Systems <u>http://www.opera.com</u>, <u>http://</u> <u>www.mozilla.com</u>. For Mac OS X, browser such as Safari, Omniweb must have an SVG plug-in installed. The plugin can be downloaded from <u>http://www.adobe.com/svg</u>. On Microsoft Windows system, a built-in web browser, IE, would also require SVG plug-in to be installed. Without the plug-in rendering all other features will still work except SVG visualization.

Abbreviations

AS-PCR, Allele Specific Polymerase Chain Reaction

ARMS, Amplification Refractory Mutation System

AS Primer, Allele Specific Primer

SNP, Single Nucleotide Polymorphism

Authors' contributions

PW and JC wrote the core engine of WASP while KC and UR contributed to the SVG visualization of the tool. CN and ST took care of all local SNP database construction. ST and AA initiated the idea of WASP. ST and PW wrote the manuscript whereas AA reviewed and tested the software.

Acknowledgements

This work is partially supported by BIOTEC grant and NRC targeted research fund. The authors wish to thank Dr. Chanin Limwongse and his members from the division of molecular genetics, Siriraj Hospital for comments and suggestions to improve this software. We also thank Suparat Chuechote for giving useful comments for improving this manuscript.

References

- Roberts RL, Kennedy MA: Rapid detection of common cytochrome p450 2d6 alleles in caucasian. Clinica Chimica Acta 2006, 366(1-2):348-51.
- 2. Bhardwaj U, Zhang YH, Lory F, McCabe LL, McCabe ERB: Molecular genetic confirmatory testing from newborn screening samples for the common African-American, Asian Indian, Southeast Asian and Chinese Thalassemia mutations. *Am J Hematol* 2005, **78(4):**249-255.
- Hasan KN, Shoji M, Tsutaya S, Kudo R, Matsuda E, Saito J, Kimura T, Yasujima M: Study of VI a vasopressin receptor gene single nucleotide polymorphisms in platelet vasopressin responsiveness. J Clin Lab Anal 2006, 20(3):87-92.
- Newton CR, Graham A, Heptinstall LE, Powell SJ, Summers C, Kalsheker N, Smith JC, Markham AF: Analysis of any point mutation in DNA. The amplification refractory mutation system (ARMS). Nucleic Acids Res 1989, 17(7):2503-16.
- Ye S, Dhillon S, Ke X, Collins AR, Day IN: An efficient procedure for genotyping single nucleotide polymorphisms. Nucleic Acids Res 2001, 29(17):E88-8.
- 6. Primo SNP, Chang BioScience [http://www.changbio_science.com/primo/primosnp.html]
- 7. SantaLucia J Jr: Physical Principles and Visual -OMP Software for Optimal PCR Design. In Method Mol Cell Biol: PCR Primer Design

Edited by: Anton Yuryev. Humana Press, Totowa, New Jersey; 2006 in press.

- 8. Little S: Amplification-Refractory Mutation System (ARMS): Analysis of Point Mutations. Curr Prot Hum Genet 1995:9.8.1-9.8.12.
- Sherry ST, Ward MH, Kholodov M, Baker J, Phan L, Smigielski EM, Sirotkin K: dbSNP: the NCBI database of genetic variation. Nucleic Acids Res 2001, 29(1):308-311.
- International HapMap Consortium: A haplotype map of the human genome. Nature 2005, 437:1299-1320.
- Hirakawa M, Tanaka T, Hashimoto Y, Kuroda M, Takagi T, Nakamura Y: JSNP: a database of common gene variations in the Japanese population. Nucleic Acids Res 2002, 30(1):158-62.
- 12. Kent J: BLAT-The BLAST-Like Alignment Tool. Genome Res 2002, 12:656-664.
- Rozen S, Skaletsky HJ: Primer3 on the WWW for general users and for biologist programmers. Methods Mol Biol 2000, 132:365-86.
- Marez D, Legrand M, Sabbagh N, Guidice JM, Spire C, Lafitte JJ, Meyer UA, Broly F: Polymorphism of the cytochrome P450 CYP2D6 gene in a European population: characterization of 48 mutations and 53 alleles, their frequencies and evolution. *Pharma*cogenetics 1997, 7(3):193-202.
- 15. Andrew W, Lilley C: SVG Unleashed. Indianapolis: SAMS 2003.

